PTS - Batchmixer

PTS-Batchmixer	Batch Size / Liters
80	50-100
100	100-250
150	250-500
200	500-1000
300	1000-2000

E With the PTS-Batchmixer it is easy to mix powders with different characteristics in a contained manner. This system is especially suited to the needs of pharmaceutical applications where active substances have to be handled without any modification of their physical properties and where all risk of contamination has to be avoided.

The system comprises a main container with an integrated central deflector. A PTS powder transfer system with two tangential inlets is installed on top of the container.

The powders are introduced automatically by the PTS-System and then circulated within the container for a predetermined period of time. A homogenous distribution of the mixture in the tank is guaranteed thanks to the deflector.

The mixing effect in the PTS body, when the two jets of powder meet, makes it possible to increase the speed and effectiveness of the mixing process significantly. Since the speed at which powders are circulated is limited, particles are not damaged. The system can operate under inert conditions, which makes it possible to mix hygroscopic, oxygen sensitive or explosive powders.

This technology is easy to integrate into pharmaceutical production lines. Powders can be transferred automatically from drums, bags or even directly from process equipment (e.g. a granulator). At the end of the mixing process the system can be emptied completely and automatically towards the next production step.

The system has no moving or rotating mechanical parts therefore, it is easy to clean automatically.

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
Special Characteristics of the PTS-Batchmixer

- A closed self-filling mixing system ensuring a high level of containment.
- Highly efficient mixing with reduced mixing times compared to traditional systems.
- High turn down ratio < 10% - 100% product load.
- Mixes diverse powders with large differences in ratio (1/10,000).
- Operates under inert conditions with low nitrogen consumption.
- The system discharges fully and can be cleaned in place (CIP).
- Small quantities of additional components can be introduced into the circulating powder without interrupting the mixing process.
- Product characteristics are not modified during mixing.
- Fixed installation which can suck powders from different types of containers (drums, Big-Bags, silos, etc.).
- Powders can be transferred over considerable distances.
- Highly reliable system with few moving parts; requires only a minimum of maintenance.

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